

OCULUS
SURGICAL

BIOM[®] Optic Set

High-definition wide field front lens
and reduction lens in single-use design



The perfect view
in every case

- Extremely wide field of view with outstanding resolution
- HD clarity under high magnification reduces the need for a contact lens
- Full-field clarity for decreased scleral depressing and panretinal laser
- Excellent depth of field for better stereopsis
- Improved view during air/fluid exchanges

The NEW BIOM® Optic Set

- Always ready, convenient blister pack design
- Minimizes risk of infection and cross contamination
- No sterilization "down time" for O.R. increased efficiency




Convenient and sterile blister package*



BIOM® HD
Disposable Lens



New reduction lens
increases surgical efficiency



BIOM® Optic Set Options

Objective lens	BIOM®	BIOM® Optic Set (disposable)
f = 175 mm	BIOM® 3 m/c, 4 m/c, 5 m/c	HD Professional for f = 175 mm
f = 200 mm	BIOM® 3 ml/cl, 4 ml/cl, 5 ml/cl	HD Professional for f = 200 mm
	BIOM® 3 m/c, 4 m/c, 5 m/c	HD Flex

* Carton contents: Single-use reduction lens, thumbscrew cover, drive belt for motorized BIOM®, BIOM® HD Disposable Lens.

CE in accordance with Medical Device Directive 93/42/EEC

WWW.OCULUSSURGICAL.COM

OCULUS Optikgeräte GmbH
Postfach • 35549 Wetzlar • GERMANY
Tel. +49 641 2005-0 • Fax +49 641 2005-299
Email: support@oculussurgical.de • www.oculussurgical.com

- OCULUS Surgical USA, info@oculussurgical.com
- OCULUS Canada, sales@oculus.ca
- OCULUS Asia, info@oculus.hk
- OCULUS Iberia, info@oculus.es
- OCULUS Turkey, info@oculus-turkey.com.tr
- OCULUS Czechia, oculus@oculus.cz



OCULUS is certified by TÜV according to
DIN EN ISO 13485

